


# Scanning Near-Field Optical Microscope (NSOM, SNOM)

## Instrument description

SNOM is a microscopic tool, which breaks the far-field light resolution limit by mapping the near-field light (evanescent waves) distribution of nanostructures. In order to achieve this, there is a very sharp optical probe (detectors/illuminators), which has an aperture of tens of nanometers. This scans the sample surface at a distance much smaller than the wavelengths of the light interacting with it. The optical resolution and structure size resolution is limited by the probe size, not by the wavelength of the incidental light (resolution  $\sim 100$  nm). This technique provides for the capabilities of the basic experimental setup – illumination by a SNOM probe set to collect interacting light in reflection or transmission; to illuminate a sample in transmission or reflection and collection the interacting light via a SNOM probe.



**Instrument:** Nanonics MultiView 4000

multi-probe independent Scanning Probe Microscope (SPM) fully integrated with upright and inverted optical microscope

**Features:**

- SPM techniques - SNOM, Atomic Force Microscopy (AFM), conductive AFM, Scanning Tunneling Microscopy (STM)
- two-independent probe scanning system, sample scanning
- combined SNOM collection, illumination with reflection and transmission modes
- probes - bent optical fibres on tuning fork according to applications
- liquid cell
- optical and acoustic hoods

## Application

	Lithographic structure	Hexacomb structure	Square of three slits structures
Topography images			
SNOM images	<p>illumination reflection</p> <p>interaction of near-field illumination with nanometer scale structure</p>	<p>collection transmission</p> <p>binary sample - transparent (centers) opaque (corners) features</p>	<p>collection reflection</p> <p>four-way evanescent waves interaction</p>

## Technical specification

### SNOM modes:

illumination transmission/reflection  
collection transmission/reflection

### AFM non-contact mode

### AFM conductive

### STM

#### Lasers:

fiber coupler + bandpass filters  
Nd:YAG  $\lambda = 532$  nm - green,  
power 20 mW, PGL-020-11-A  
HeNe  $\lambda = 632.8$  nm - red,  
power 10 mW, JDSU 1135/P

#### Detectors:

Avalanche PhotoDiode (APD)  
- SPCM-AQR-14 Perkin Elmer  
-  $\lambda \in <400, 1100>$  nm  
- photoncounting mode  
PhotoMultiplier Tube (PMT)  
- MP942 Perkin Elmer  
-  $\lambda \in <165, 650>$  nm  
- photoncounting mode


#### Optical microscopes:

upright - Olympus BXF  
inverted - Olympus BXF  
objectives - 50x NA 0.45, 10x  
confocal input/output module - pinhole 50  $\mu$ m  
Attenuated Total Reflection module

#### Manipulators:

xyz scanner 1: 40  $\mu$ m x 40  $\mu$ m x 30  $\mu$ m  
xyz scanner 2: 40  $\mu$ m x 40  $\mu$ m x 30  $\mu$ m  
xyz sample scanner: 80  $\mu$ m x 80  $\mu$ m x 30  $\mu$ m  
xyz micromanipulators: 5 mm x 5 mm x 10 mm

fiber probes optimized for selected measurement



## Contact

**Core Facility:** Nanofabrication and Nanocharacterization

**Section:** Probe Microscopy & Nanomanipulation

**Contact person:** David Škoda

[core.facility@ceitec.vutbr.cz](mailto:core.facility@ceitec.vutbr.cz)

**Detailed information:** [www.ceitec.eu](http://www.ceitec.eu)

**Instrument location:**

Masaryk University  
Faculty of Science  
Department of Condensed Matter Physics  
Kotlářská 267/2, 611 37 Brno  
Czech Republic